

D95 SET

WATER IN PETROLEUM PRODUCTS AND BITUMINOUS MATERIALS BY DISTILLATION

SPECIFICATIONS

Standards

ASTM D 95, IP 74, ISO 3733 and related methods.

Key features

- Electric mantle and energy regulator

Use

The amount of water determined by this test method may be used to correct the volume involved in the custody transfer of petroleum products and bituminous materials.



ORDERING INFORMATION

Scope of delivery :

21454 Instrument is delivered with:

- 1 heating mantle PN 26402
- 1 round bottom flask 500 ml with 24/29 ground joint PN 12852
- 1 Dean Stark (receiver) round bottom 5 ml, with ground joint 24/29 PN 21455
- 1 graduated Liebig Condenser 1/10, 400 mm effective length with 24/29 ground joint PN 13142
- 1 stand PN 16324
- 1 boss PN 16329
- 1 clamp PN 16252
- 1 lift holder PN 11048
- 1 silicone tube PN 14094

Site requirements:

Power supply: 220V – 50 Hz



CONTACT :
sales@normalab.com (Export)

Normalab France SAS
76190 Valliquerville - France
Tel.: +33 232.700.100
Fax.: +33 232.704.732

DISTRIBUTED BY